

## pH FIELD & LAB ELECTRODES

Clear Epoxy-Bodied, Gel-Filled Combination ALphA® Electrodes

PHE-4201



PHE-4232 semi-micro.

PHE-4201 general-purpose.

PHE-4202

double junction.



PHE-4272 flat surface double junction.

✓ Five Electrode Styles

connector standard.

- Covers Entire 0 to 14 pH Range
- Clear Epoxy Body Construction
- ✓ Fast Responding: 95% in Less than 1 Second
- ✓ Comes with OMEGA **Quality Certificate**

OMEGA Engineering is pleased to offer a complete line of laboratory pH electrodes. The PHE-4200 Series consists of gel-filled, combination pH electrodes with a unique clear epoxy body design. With five electrode styles to choose from, you're sure to find one suitable for your application.

OMEGA® PHE-4200 Series electrodes offer reliable quality at a reasonable cost. The electrodes measure the entire pH range from 0 to 14 pH units at 0 to 100°C (32 to 212°F) and feature an Ag/AgCI reference and polypropylene liquid junction. Each electrode comes in a soaker storage bottle to keep the electrode moist and ready for use. Shelf life is guaranteed for one year, and the provided OMEGA quality certificate assures that these electrodes meet the most stringent quality control requirements.


## Certificate of Quality

PHE-4222

flask size.

This electrode has been manufactured to the most stringent quality control standards. It has been shipped in a special soaker/storage bottle to ensure its quality.

Your electrode has been 100% inspected prior to shipment and meets or exceeds all factory performance specifications.



To Order				
Model No.	Description	Length mm (inch)	Dia. mm (inch)	R (25°C) ΜΩ
PHE-4201	General purpose	150 (6)	12 (0.47)	50 to 80
PHE-4202	Double junction	150 (6)	12 (0.47)	50 to 80
PHE-4222	Flask size	300 (12)	9.5 (0.37)	50 to 80
PHE-4272	Flat surface double junction	150 (6)	12 (0.47)	50 to 80
PHE-4232	Semi-micro	150 (6)	6 (0.24)	50 to 80
PHA-D1	500 ml (1 pint) bottle deionized water			

All electrodes come complete with 750 mm (2.5') of cable, BNC connector and electrode care instruction sheet.

Ordering Examples: PHE-4201, general purpose electrode. PHE-4202, double junction electrode.